

Disposable COOLVEST

Made in the U.S.A

Instructions for Use:

- 1) Fill cooler with ice and 1 gallon of water
 - For use with Cooling Pad skip to #4
- 2) Put on Disposable COOLVEST®
- 3) Adjust Velcro as necessary.
- 4) Attach hoses to Disposable COOLVEST® or Cooling Pad. Audible click means secured hoses.
- 5) Press ON/OFF switch to ON

NOTES:

- Water hose connections are quick release dry-break connectors
- Upon initial use, allow 1-2 minutes for the Disposable COOLVEST® to fill with water
- Use hose flow controller to control temperature of Disposable COOLVEST®

Troubleshooting:

- 1) **Pump Motor Not Running:** Check electrical connections.
- 2) **Water not flowing:** Reconnect hoses at cooler and Disposable COOLVEST®
Add water if needed

Maintenance:

COOLER:

- Add 2 ounces of Maintenance Additive (MA-16) to ice water with each use.
- Dry inside of cooler with a soft cloth after each use.

HOSES:

- Replace O-rings yearly or as needed.
- Use a soft cloth with mild soap and water to clean hoses.